



Product description

Order number	65.386.007
Product description	Colour-coded inserts, for CryoPure tubes, violet, 100 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	CryoPure tubes

Material & colours

Product material	Polystyrene (PS)
Colour of product	violet

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	3000
Piece(s) / pallet	972000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0002 kg
Weight of case	0.68 kg
EAN of inner box	4038917106184
EAN case	4038917077729

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

