



### **Product description**

Order number	72.379.004
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, yellow, external thread, nominal volume: 1.8 ml, Cryo Performance Tested, 50 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	external thread
Volume of work	1,8 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**CryoPure tubes, 2 ml, QuickSeal screw cap,  
yellow**

**Page 2**

**Size**

Sample volume	1.8 ml
Max. volume	2 ml
Diameter	12.65 mm
Length of product	46.7 mm
Length excluding cap	39.8 mm

**Material & colours**

Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.002 kg
Weight of case	5.17 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-04-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**CryoPure tubes, 2 ml, QuickSeal screw cap,  
yellow**

**Page 3**

EAN of inner box                      4038917106047  
EAN case                                 4038917080576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-04-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**