



**Product description**

Order number	72.380.007
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, PP, violet, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	Internal thread
Design	with silicon O-ring
Volume of work	1,6 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Sample volume	1.6 ml
Max. volume	2 ml
Diameter	12.65 mm
Length of product	47.7 mm
Length excluding cap	37.8 mm

**Material & colours**

Closure material	Polypropylene (PP)
Colour of cap	violet
O-ring material	silicone
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0019 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**CryoPure tubes, 2 ml, QuickSeal screw cap, violet**

Weight of case	4.99 kg
EAN of inner box	4038917106108
EAN case	4038917080644

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

