



Product Specification
**CryoPure tubes, 5 ml, QuickSeal screw cap,
yellow**

Page 1



Product description

Order number	72.383.004
Product description	CryoPure tubes, 5 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, yellow, external thread, nominal volume: 4.5 ml, Cryo Performance Tested, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	external thread
Volume of work	4,5 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**CryoPure tubes, 5 ml, QuickSeal screw cap,
yellow**

Page 2

Size

Sample volume	4.5 ml
Max. volume	5 ml
Diameter	12.65 mm
Length of product	90.7 mm
Length excluding cap	83.8 mm

Material & colours

Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	52000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0038 kg
Weight of case	5.05 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**CryoPure tubes, 5 ml, QuickSeal screw cap,
yellow**

Page 3

EAN of inner box	4038917106122
EAN case	4038917080675

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT