



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 93.874.410  |
| Product description | Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, yellow, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 10 x 10, for 100 collection tubes, for CryoPure tubes 1.2?2.0 ml internal thread, 5 piece(s)/bag |

**Product characteristics**

|                  |   |
|------------------|---|
| Graduation       | no  |
| Cap              | slip-on lid with ventilation function     |
| Closure type     | slip-on lid                               |
| Design           | with numerical coding at each aperture    |
| Suitable for     | CryoPure tubes 1.2?2.0 ml internal thread |
| Storage capacity | 100                                       |
| Format           | 10 x 10                                   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Size**

|                   |        |
|-------------------|--------|
| Width of product  | 132 mm |
| Height of product | 53 mm  |
| Length of product | 132 mm |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polycarbonate (PC) |
| Colour of product | yellow             |
| Colour of cap     | transparent        |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5             |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 5             |
| Piece (s) / outer case        | 20            |
| Piece(s) / pallet             | 720           |
| Depth of case                 | 135 mm        |
| Width of case                 | 270 mm        |
| Height of case                | 540 mm        |
| Case volume                   | 0.0197 cbm    |
| Weight of product             | 0.1593 kg     |
| Weight of case                | 3.71 kg       |
| EAN of inner box              | 4038917122917 |
| EAN case                      | 4038917415163 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

