




**Product Specification**  
**Tissue culture dish, (ØxH): 60 x 15 mm, surface:**  
**Suspension**

**Page 1**



**Product description**

Order number	83.3901.500
Product description	Tissue culture dish, (ØxH): 60 x 15 mm, material: PS, surface: Suspension, for suspension cells, coding colour: green, TC Tested, 10 piece(s)/bag

**Product characteristics**

Closure type	slip-on lid
Cell type	suspension cells
Base shape	SUREGrip Design
Growth area	22.1 cm <sup>2</sup>
Volume of work	5 ml
Vessel type	Dish
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tissue culture dish, (ØxH): 60 x 15 mm, surface:**  
**Suspension**

**Page 2**

**Size**

Diameter	60 mm
Height of product	15 mm

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	332 mm
Width of case	332 mm
Height of case	380 mm
Case volume	0.0419 cbm
Weight of product	0.0079 kg
Weight of case	5.02 kg
EAN of inner box	4038917338134
EAN case	4038917338127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

