



**Product Specification**  
**Cell culture flask, T-25, surface: Standard,  
2-position screw cap**

**Page 1**



**Product description**

Order number	83.3910
Product description	Cell culture flask, T-25, material: PS, surface: Standard, for Adherent cells, coding colour: red, 2-position screw cap, TC Tested, 10 piece(s)/bag

**Product characteristics**

Cap	2-position screw cap
Cell type	Adherent cells
Growth area	25 cm <sup>2</sup>
Volume of work	7 ml - 12,5 ml
Vessel type	flask
Coding colour	red

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cell culture flask, T-25, surface: Standard,**  
**2-position screw cap**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	300
Piece(s) / pallet	6000
Depth of case	602 mm
Width of case	408 mm
Height of case	340 mm
Case volume	0.0835 cbm
Weight of product	0.0142 kg
Weight of case	6.61 kg
EAN of inner box	4038917337595
EAN case	4038917337588

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**