



### **Product description**

Order number	83.3911.300
Product description	Cell culture flask, T-75, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, 2-position screw cap, TC Tested, 5 piece(s)/bag

### **Product characteristics**

Cap	2-position screw cap
Cell type	challenging adherent cells
Growth area	75 cm <sup>2</sup>
Volume of work	21 ml - 55 ml
Vessel type	flask
Coding colour	yellow

### **Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1600
Depth of case	445 mm
Width of case	402 mm
Height of case	435 mm
Case volume	0.0778 cbm
Weight of product	0.0461 kg
Weight of case	6.6 kg
EAN of inner box	4038917337472
EAN case	4038917366069

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

