



Product description

Order number	83.3912.500
Product description	Cell culture flask, T-175, material: PS, surface: Suspension, for suspension cells, coding colour: green, 2-position screw cap, TC Tested, 5 piece(s)/bag

Product characteristics

Cap	2-position screw cap
Cell type	suspension cells
Growth area	175 cm ²
Volume of work	50 ml - 125 ml
Vessel type	flask
Coding colour	green

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	40
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	40
Piece(s) / pallet	960
Depth of case	522 mm
Width of case	267 mm
Height of case	450 mm
Case volume	0.0627 cbm
Weight of product	0.106 kg
Weight of case	5.63 kg
EAN of inner box	4038917337892
EAN case	4038917337885

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

