



**Product Specification**  
**Cell culture flask, T-175, surface: Standard,  
2-position screw cap**

**Page 1**



**Product description**

Order number	83.3912
Product description	Cell culture flask, T-175, material: PS, surface: Standard, for Adherent cells, coding colour: red, 2-position screw cap, TC Tested, 5 piece(s)/bag

**Product characteristics**

Cap	2-position screw cap
Cell type	Adherent cells
Growth area	175 cm <sup>2</sup>
Volume of work	50 ml - 125 ml
Vessel type	flask
Coding colour	red

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cell culture flask, T-175, surface: Standard,**  
**2-position screw cap**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	40
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	40
Piece(s) / pallet	960
Depth of case	522 mm
Width of case	267 mm
Height of case	450 mm
Case volume	0.0627 cbm
Weight of product	0.106 kg
Weight of case	5.63 kg
EAN of inner box	4038917337717
EAN case	4038917337700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-24

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**