



### **Product description**

Order number	83.3924.005
Product description	Cell culture plate, 96 well, material: PS, surface: Standard, for Adherent cells, coding colour: red, flat base, TC Tested, 5 piece(s)/bag

### **Product characteristics**

Label/ Print	with plastic label
Closure type	slip-on lid
Cell type	Adherent cells
Base shape	flat base
Number of wells	96
Growth area	0.29 cm <sup>2</sup>
Volume of work	0,2 ml
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Max. volume	0.39 ml
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	2400
Depth of case	567 mm
Width of case	372 mm
Height of case	280 mm
Case volume	0.0591 cbm
Weight of product	0.0568 kg
Weight of case	8.17 kg
EAN of inner box	4038917360258
EAN case	4038917360241

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

