



### **Product description**

Order number	83.3945.999
Product description	Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### **Product characteristics**

Suitable for	15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes
--------------	-----------------------------------------------

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm**  
**and 12 x 75 mm tubes**

**Page 2**

**Packaging**

Minimum order qty.	25
Type of smallest subpackaging	blister
Piece(s) / inner box	25
Piece (s) / outer case	200
Piece(s) / pallet	8400
Depth of case	229 mm
Width of case	330 mm
Height of case	584 mm
Case volume	0.0441 cbm
Weight of product	0.0053 kg
Weight of case	2 kg
EAN of inner box	4038917415323
EAN case	4038917415316

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**