




**Product Specification**  
**x-well cell culture chamber, 4 wells, on lumox® slide, removable frame**

**Page 1**



**Product description**

Order number	94.6150.401
Product description	x-well cell culture chamber, 4 wells, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

**Product characteristics**

Closure type	slip-on lid
Design	on lumox® slide
Cell type	Adherent cells
Number of wells	4
Growth area	1.9 cm <sup>2</sup>
Volume of work	1 ml
Removable frame	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**x-well cell culture chamber, 4 wells, on lumox®**  
**slide, removable frame**

**Page 2**

**Size**

Width of product	24 mm
Length of product	76 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	6
Type of smallest subpackaging	blister
Piece(s) / inner box	6
Piece (s) / outer case	96
Piece(s) / pallet	10368
Depth of box	127 mm
Width of box	19 mm
Height of box	258 mm
Depth of case	275 mm
Width of case	270 mm
Height of case	185 mm
Case volume	0.0137 cbm
Weight of product	0.0472 kg
Weight of case	5.41 kg
EAN of inner box	4038917255622
EAN case	4038917255615

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**