



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 83.1826.102   |
| Product description | Syringe filter, Filtropur S plus, pre-filter: GF, membrane: CA, filtration surface: 6.2 cm², pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

**Product characteristics**

|                        |                        |
|------------------------|------------------------|
| Example of use         | for sterile filtration |
| Filter                 | yes                    |
| Surface area of filter | 6.2 cm²                |
| Pore size              | 0,2                    |
| Connection             | Luer lock              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Syringe filter, Filtropur S plus, CA, pore size: 0.2 µm, for sterile filtration**

**Page 2**

**Size**

Diameter 33 mm

**Material & colours**

Product material Methacrylate-butadiene-styrene (MBS)  
Colour of product light blue  
Membrane material Cellulose acetate (CA)

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic  
Sterilisation Gamma irradiation  
Pyrogen-free yes  
Non-cytotoxic yes  
Batched yes

**Packaging**

Minimum order qty. 50  
Type of smallest subpackaging blister  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 12000  
Depth of case 375 mm  
Width of case 285 mm  
Height of case 530 mm  
Case volume 0.0566 cbm  
Weight of product 0.0064 kg  
Weight of case 4.2 kg  
EAN of inner box 4038917133920  
EAN case 4038917025263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

