



Product description

| | |
|---------------------|---|
| Order number | 83.3942.100 |
| Product description | Filtropur BT 100, Bottle top filter, for screw cap 45 mm, 1,000 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |

Product characteristics

| | |
|------------------------|----------------------|
| Graduation | yes |
| Example of use | for clear filtration |
| Suitable for | screw cap 45 mm |
| Filter | yes |
| Surface area of filter | 63.6 cm² |
| Pore size | 0,45 |
| Connection | Step taper |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------------|
| Sample volume | 1.000.00 ml |
| Filtration volume | 1000 ml |
| Height of product | 15.27 cm |

Material & colours

| | |
|-------------------|------------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Membrane material | Polyethersulfone (PES) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 24 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 24 |
| Piece(s) / pallet | 288 |
| Depth of case | 572 mm |
| Width of case | 343 mm |
| Height of case | 439 mm |
| Case volume | 0.0861 cbm |
| Weight of product | 0.1305 kg |
| Weight of case | 5.08 kg |
| EAN of inner box | 4038917412322 |
| EAN case | 4038917412315 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

