




Product Specification
**Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm
and 12 x 75 mm tubes**

Page 1



Product description

Order number	83.3945.999
Product description	Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Suitable for	15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes
--------------	---

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm
and 12 x 75 mm tubes

Page 2

Packaging

Minimum order qty.	25
Type of smallest subpackaging	blister
Piece(s) / inner box	25
Piece (s) / outer case	200
Piece(s) / pallet	8400
Depth of case	229 mm
Width of case	330 mm
Height of case	584 mm
Case volume	0.0441 cbm
Weight of product	0.0053 kg
Weight of case	2 kg
EAN of inner box	4038917415323
EAN case	4038917415316

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT