




**Product Specification**  
**lumox® dish 35, Tissue culture dish, with foil**  
**base, Ø: 35 mm, suspension cells**

**Page 1**



**Product description**

Order number	94.6077.333
Product description	lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

**Product characteristics**

Closure type	slip-on lid
Design	with foil base
Cell type	suspension cells
Base shape	flat base
Growth area	6.35 cm <sup>2</sup>
Volume of work	2,5 ml
Vessel type	Dish
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**lumox® dish 35, Tissue culture dish, with foil  
base, Ø: 35 mm, suspension cells**

**Page 2**

**Size**

Diameter	35 mm
Height of product	6 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	15000
Depth of case	348 mm
Width of case	253 mm
Height of case	265 mm
Case volume	0.0233 cbm
Weight of product	0.0037 kg
Weight of case	1.91 kg
EAN of inner box	4038917255141
EAN case	4038917255134

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**