




Product Specification
lumox® multiwell, Cell culture plate, with foil
base, 24 well, 20 piece(s)

Page 1



Product description

Order number	94.6110.024
Product description	lumox® multiwell, Cell culture plate, with foil base, 24 well, frame: black, frame: PS, for Adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Closure type	slip-on lid
Design	with foil base
Cell type	Adherent cells
Base shape	flat base
Number of wells	24
Growth area	1.9 cm ²
Volume of work	500 µl - 1.500 µl
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**lumox® multiwell, Cell culture plate, with foil
base, 24 well, 20 piece(s)**

Page 2

Size

Width of product	85.48 mm
Height of product	23.4 mm
Length of product	127.76 mm

Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)
Colour of frame	black

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1920
Depth of case	333 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0151 cbm
Weight of product	0.065 kg
Weight of case	1.67 kg
EAN of inner box	4038917255448
EAN case	4038917255431

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

