



### **Product description**

Order number	94.6001.054
Product description	miniPERM® , Rack, for miniPERM® bioreactor, not autoclavable, 4 piece(s)/case

### **Product characteristics**

Suitable for	miniPERM® bioreactor
--------------	----------------------

### **Material & colours**

Product material	Acrylnitrile butadiene styrene (ABS)
Colour of product	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	4
Type of smallest subpackaging	case
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	672
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.085 kg
Weight of case	3.7 kg
EAN of inner box	4038917254786
EAN case	4038917254779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

