



Product description

| | |
|---------------------|---|
| Order number | 94.6001.055 |
| Product description | miniPERM® , classic production module, Tissue culture compartment, for protein, virus and biomass production, for suspension cells, volume of product module: 35 ml, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |

Product characteristics

| | |
|--------------|---|
| Closure type | screw cap |
| Design | for suspension cells |
| Application | for protein, virus and biomass production |
| Cell type | suspension cells |

Size

| | |
|--------------------------|--------|
| Production module volume | 35 ml |
| Diameter | 102 mm |
| Height of product | 37 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|----------------------------|------------------------|
| Product material | Polycarbonate (PC) |
| Closure material | Polypropylene (PP) |
| Dialysis membrane material | Cellulose acetate (CA) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 4 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 4 |
| Piece (s) / outer case | 24 |
| Piece(s) / pallet | 480 |
| Depth of case | 498 mm |
| Width of case | 364 mm |
| Height of case | 310 mm |
| Case volume | 0.0562 cbm |
| Weight of product | 0.1414 kg |
| Weight of case | 5.13 kg |
| EAN of inner box | 4038917385619 |
| EAN case | 4038917385596 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

