



## Product description

Order number	94.6001.094
Product description	miniPERM® , Start-up support kit, accessories for inoculation, sample collection and harvesting, for miniPERM® bioreactor, not autoclavable, set components: 8 x disposable syringe 50 ml, 20 x disposable syringe 2 ml, 8 x filling tube, 20 x needle, 6 x septum cap, 1 x miniPERM® stand, 1 piece(s)/bag

## Product characteristics

Application	accessories for inoculation, sample collection and harvesting
Suitable for	miniPERM® bioreactor

## Purity & certification

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	96
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.481 kg
Weight of case	0.93 kg
EAN of inner box	4038917254939
EAN case	4038917254922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

