




**Product Specification**  
**x-well cell culture chamber, 1 well, on lumox®**  
**slide, removable frame**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 94.6150.101   |
| Product description | x-well cell culture chamber, 1 well, on lumox® slide,<br>removable frame, sterile, non-pyrogenic/endotoxin-free,<br>non-cytotoxic, 6 piece(s)/blister |

**Product characteristics**

|                 |                   |
|-----------------|-------------------|
| Closure type    | slip-on lid       |
| Design          | on lumox® slide   |
| Cell type       | Adherent cells    |
| Number of wells | 1                 |
| Growth area     | 9 cm <sup>2</sup> |
| Volume of work  | 4 ml              |
| Removable frame | yes               |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**x-well cell culture chamber, 1 well, on lumox®**  
**slide, removable frame**

**Page 2**

**Size**

|                   |       |
|-------------------|-------|
| Width of product  | 24 mm |
| Length of product | 76 mm |

**Material & colours**

|                  |                  |
|------------------|------------------|
| Product material | lumox®           |
| Frame material   | Polystyrene (PS) |

**Purity & certification**

|                  |  |
|------------------|--|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gamma irradiation                                    |
| Pyrogen-free     | yes  |
| Batched          | yes  |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 6             |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 6             |
| Piece (s) / outer case        | 96            |
| Piece(s) / pallet             | 10368         |
| Depth of box                  | 127 mm        |
| Width of box                  | 19 mm         |
| Height of box                 | 258 mm        |
| Depth of case                 | 275 mm        |
| Width of case                 | 270 mm        |
| Height of case                | 185 mm        |
| Case volume                   | 0.0137 cbm    |
| Weight of product             | 0.0461 kg     |
| Weight of case                | 5.3 kg        |
| EAN of inner box              | 4038917255561 |
| EAN case                      | 4038917255554 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

