



## Product description

Order number 90.183.200  
Product description adapter, for centrifuge SC 2700

## Product characteristics

Suitable for centrifuge SC 2700

## Size

Diameter 17 mm  
Width of product 17 mm  
Height of product 12.5 mm  
Length of product 17 mm

## Material & colours

Product material Polypropylene (PP)  
Colour of product grey

## Purity & certification

Product category no medical device | IVD  
Batched no

## Packaging

Minimum order qty. 2  
Type of smallest subpackaging bag  
Piece(s) / inner box 2  
Piece (s) / outer case 300  
Piece(s) / pallet 97200  
Depth of case 268 mm  
Width of case 191 mm  
Height of case 100 mm  
Case volume 0.0051 cbm  
Weight of product 0.002 kg  
Weight of case 1 kg  
EAN of inner box 4038917909679  
EAN case 4038917909686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

