



Product description

Order number 90.2620.304
Product description length: 1.5 m, transparent, 1 piece(s)/bag

Product characteristics

Suitable for IVARO

Size

Diameter 0.8 mm
Length of product 1.5 m

Material & colours

Product material silicone
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging bag
Piece(s) / inner box 1
Piece (s) / outer case 10
Case volume 0.0062 cbm
Weight of product 0.0005 kg
EAN of inner box 4038917584111
EAN case 4038917584104

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

