



Product description

| | |
|---------------------|---|
| Order number | 92.2620.304 |
| Product description | length: 1.5 m, material: silicone, transparent, for IVARO, 1 piece(s)/bag |

Product characteristics

| | |
|--------------|-------|
| Suitable for | IVARO |
|--------------|-------|

Size

| | |
|-------------------|--------|
| Diameter | 0.8 mm |
| Length of product | 1.5 m |

Material & colours

| | |
|-------------------|-------------|
| Product material | silicone |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 10 |
| Depth of case | 280 mm |
| Width of case | 220 mm |
| Height of case | 100 mm |
| Case volume | 0.0062 cbm |
| Weight of product | 0.0005 kg |
| EAN of inner box | 4038917584111 |
| EAN case | 4038917584104 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

