




**Product Specification**  
**Deep Well MegaBlock® , 0.5 ml, Free from**  
**DNA/DNase/RNase, free from**  
**pyrogens/endotoxins, PP**

**Page 1**



**Product description**

Order number	82.1969.002
Product description	Deep Well MegaBlock® , 96 well, 0.5 ml, round base, Free from DNA/DNase/RNase, free from pyrogens/endotoxins, material: PP, SLAS standard, Round cavities, 7 piece(s)/bag

**Product characteristics**

Closure type	film
Design	SLAS standard
Application	collection and storage of diagnostic specimens
Base shape	round base
Number of wells	96
Volume of work	1 µl - 0,5 ml
Centrifugation max (RCF)	2000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Deep Well MegaBlock® , 0.5 ml, Free from**  
**DNA/DNase/RNase, free from**  
**pyrogens/endotoxins, PP**

**Page 2**

**Size**

Sample volume	0.5 ml
Diameter	7.2 mm
Width of product	85.1 mm
Height of product	28.5 mm
Length of product	127.3 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	Free from DNA/DNase/RNase, free from pyrogens/endotoxins
Batched	yes

**Packaging**

Minimum order qty.	56
Type of smallest subpackaging	bag
Piece(s) / inner box	7
Piece (s) / outer case	56
Piece(s) / pallet	4480
Case volume	0.0181 cbm
Weight of product	0.057 kg
Weight of case	3.74 kg
EAN of inner box	4038917346320
EAN case	4038917346313

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**