



**Product Specification**  
**Deep Well plate, 2.2 ml, free from**  
**DNA/DNase/RNase/PCR inhibitors, PP**

**Page 1**



**Product description**

Order number	82.1972
Product description	Deep Well plate, 96 well, 2.2 ml, compatible with KingFisher? Flex/Duo Prime/Presto/Apex, Bio Sprint 96, Chemagic? Prime?, STARlet, MultiMACS M96 Separator, conical wedge, free from DNA/DNase/RNase/PCR inhibitors, material: PP, Edged cavities, 4 piece(s)/bag

**Product characteristics**

Closure type	film
Application	collection and storage of diagnostic specimens
Suitable for	KingFisher? Flex/Duo Prime/Presto/Apex, Bio Sprint 96, Chemagic? Prime?, STARlet, MultiMACS M96 Separator
Base shape	conical wedge
Number of wells	96
Volume of work	1 µl - 2,2 ml
Centrifugation max (RCF)	4500 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Deep Well plate, 2.2 ml, free from  
DNA/DNase/RNase/PCR inhibitors, PP**

**Page 2**

**Size**

Sample volume	2.2 ml
Max. volume	2.2 ml
Diameter	7.52 mm
Width of product	85 mm
Height of product	41.85 mm
Length of product	127.76 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase/PCR inhibitors
Batched	yes

**Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	2304
Depth of case	283 mm
Width of case	268 mm
Height of case	210 mm
Case volume	0.0159 cbm
Weight of product	0.0848 kg
Weight of case	2.28 kg
EAN of inner box	4038917952477
EAN case	4038917952484

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**