




Product Specification

**Closure mat, for Deep Well MegaBlock® 0.5 ml
(82.1969.002) & 1.2 ml (82.1970.002 &
82.1971.002), (LxW): 79 x 121 mm, TPE,
transparent**



Product description

Order number	95.1990.002
Product description	Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, material: TPE, transparent, pierceable, Free from DNA/DNase/RNase, free from pyrogens/endotoxins, 10 piece(s)/bag

Product characteristics

Suitable for	Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002)
--------------	--

Size

Width of product	121 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



