



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 86.1198  |
| Product description | Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag |

### **Product characteristics**

|                        |               |
|------------------------|---------------|
| Label/ Print           | without print |
| Liquid-tight           | yes           |
| Autoclavable (max. °C) | 134 °C        |

### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 300 mm |
| Length of product | 500 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 7 l, (LxW): 500 x 300 mm, PP,**  
**transparent, without print**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                   |                         |
|-------------------|-------------------------|
| Product category  | no medical device   IVD |
| Latex-free        | yes                     |
| Safety data sheet | no                      |
| Batched           | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 48000         |
| Depth of case                 | 348 mm        |
| Width of case                 | 278 mm        |
| Height of case                | 140 mm        |
| Case volume                   | 0.0135 cbm    |
| Weight of product             | 0.0135 kg     |
| Weight of case                | 7 kg          |
| EAN of inner box              | 4038917085724 |
| EAN case                      | 4038917004350 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

