



**Product Specification**  
**Disposal bags, 40 l, (LxW): 780 x 600 mm, PP,**  
**yellow, with print**

**Page 1**



**Product description**

Order number	86.1203.103
Product description	Disposal bags, volume: 40 l, (LxW): 780 x 600 mm, material: PP, yellow, film thickness: 50 µm, autoclavable up to 134 °C, with print, label/print: blue, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	blue
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

**Size**

Width of product	600 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 40 l, (LxW): 780 x 600 mm, PP,**  
**yellow, with print**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	7000
Depth of case	623 mm
Width of case	288 mm
Height of case	115 mm
Case volume	0.0206 cbm
Weight of product	0.0415 kg
Weight of case	10.7 kg
EAN of inner box	4038917126649
EAN case	4038917020732

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**