



## Product description

Order number	86.1206.103
Product description	Disposal bags, volume: 80 l, (LxW): 1,120 x 700 mm, material: PP, yellow, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: blue, 50 piece(s)/bag

## Product characteristics

Label/ Print	with print
Colour of print/label	blue
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

## Size

Width of product	700 mm
Length of product	1.120.00 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	150
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	150
Piece(s) / pallet	3600
Depth of case	378 mm
Width of case	122 mm
Height of case	810 mm
Case volume	0.0374 cbm
Weight of product	0.072 kg
Weight of case	11.38 kg
EAN of inner box	4038917160452
EAN case	4038917160445

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

