



### **Product description**

Order number	82.1581.001
Product description	Micro test plate, 96 well, slip-on lid, flat base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

### **Product characteristics**

Label/ Print	with plastic label
Closure type	slip-on lid
Base shape	flat base
Number of wells	96
Vessel type	Plate

### **Size**

Max. volume	390 µl
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro test plate, 96 well, slip-on lid, flat base, PS,  
transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	3200
Depth of case	332 mm
Width of case	297 mm
Height of case	215 mm
Case volume	0.0212 cbm
Weight of product	0.0568 kg
Weight of case	3.54 kg
EAN of inner box	4038917131056
EAN case	4038917347952

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**