



Product description

Order number	82.1581.110
Product description	ELISA plate, 96 well, flat base, without lid, material: PS, white, Medium Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Base shape	flat base
Number of wells	96
Vessel type	Plate

Size

Max. volume	390 µl
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	white

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	50
Piece(s) / pallet	2800
Depth of case	279 mm
Width of case	232 mm
Height of case	365 mm
Case volume	0.024 cbm
Weight of product	0.0448 kg
Weight of case	3.06 kg
EAN of inner box	4038917365307
EAN case	4038917365291

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

