



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 82.1582.001  |
| Product description | Micro test plate, 96 well, slip-on lid, round base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister |

### **Product characteristics**

|                 |                    |
|-----------------|--------------------|
| Label/ Print    | with plastic label |
| Closure type    | slip-on lid        |
| Base shape      | round base         |
| Number of wells | 96                 |
| Vessel type     | Plate              |

### **Size**

|                   |          |
|-------------------|----------|
| Max. volume       | 310 µl   |
| Width of product  | 85.2 mm  |
| Height of product | 16.55 mm |
| Length of product | 127.8 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



## Material & colours

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

## Purity & certification

|                  |  |
|------------------|--|
| Product category | no medical device   IVD  |
| Purity standard  | free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile |
| Sterilisation    | Gamma irradiation  |
| Batched          | yes  |

## Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 50            |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 50            |
| Piece(s) / pallet             | 3200          |
| Depth of case                 | 332 mm        |
| Width of case                 | 297 mm        |
| Height of case                | 215 mm        |
| Case volume                   | 0.0212 cbm    |
| Weight of product             | 0.0498 kg     |
| Weight of case                | 3.19 kg       |
| EAN of inner box              | 4038917131063 |
| EAN case                      | 4038917347969 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.