



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 82.1586  |
| Product description | Adhesive film, for ELISA and microtest plates, material: acetate, transparent, self-adhesive, 100 piece(s)/bag |

### **Product characteristics**

|              |                            |
|--------------|----------------------------|
| Suitable for | ELISA and microtest plates |
|--------------|----------------------------|

### **Material & colours**

|                   |             |
|-------------------|-------------|
| Product material  | acetate     |
| Colour of product | transparent |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 100           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 30000         |
| Depth of case                 | 95 mm         |
| Width of case                 | 239 mm        |
| Height of case                | 225 mm        |
| Case volume                   | 0.0051 cbm    |
| Weight of product             | 0.0028 kg     |
| Weight of case                | 3 kg          |
| EAN of inner box              | 4038917087148 |
| EAN case                      | 4038917005609 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

