



**Product description**

Order number	71.9923.211
Product description	Weigh tray, 5 ml, (LxW): 35 x 35 mm, material: PVC, white, 200 piece(s)/bag

**Size**

Sample volume	5 ml
Max. volume	8 ml
Width of product	35 mm
Length of product	35 mm

**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	white

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	308000
Depth of case	271 mm
Width of case	158 mm
Height of case	224 mm
Case volume	0.0096 cbm
Weight of product	0.0007 kg
Weight of case	1.6 kg
EAN of inner box	4038917153652
EAN case	4038917038232

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

