



Product description

Order number	81.970
Product description	Stirring rod, length: 120 mm, material: PP, white, 1,000 piece(s)/bag

Size

Diameter	3 mm
Length of product	120 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	white

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Stirring rod, length: 120 mm, PP, white

Page 2

Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	400000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0006 kg
Weight of case	6.23 kg
EAN of inner box	4038917269186
EAN case	4038917003506

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT