



Product Specification
Stomacher bag, (LxW): 300 x 170 mm, LD-PE,
transparent

Page 1



Product description

| | |
|---------------------|---|
| Order number | 86.9924.531 |
| Product description | Stomacher bag, (LxW): 300 x 170 mm, material: LD-PE, transparent, Food testing laboratories |

Product characteristics

| | |
|----------------|---|
| Application | Food testing laboratories |
| Example of use | Prolonged massage in a Stomacher® machine |

Size

| | |
|-------------------|--------|
| Width of product | 170 mm |
| Length of product | 300 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Stomacher bag, (LxW): 300 x 170 mm, LD-PE,
transparent**

Page 2

Purity & certification

| | |
|-------------------|-----|
| Latex-free | yes |
| Pyrogen-free | no |
| Safety data sheet | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 80000 |
| Depth of case | 458 mm |
| Width of case | 238 mm |
| Height of case | 373 mm |
| Case volume | 0.0408 cbm |
| Weight of case | 9.3 kg |
| EAN of inner box | 4038917155380 |
| EAN case | 4038917398442 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT