



Product description

Order number	91.2620.020
Product description	Thermo transfer foil, black, 1 piece(s)/roll

Material & colours

Colour of product	black
-------------------	-------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	2
Type of smallest subpackaging	roll
Piece(s) / inner box	1
Piece (s) / outer case	2
Piece(s) / pallet	648
Depth of case	260 mm
Width of case	183 mm
Height of case	85 mm
Case volume	0.004 cbm
Weight of product	0.209 kg
Weight of case	0.52 kg
EAN of inner box	4038917899703
EAN case	4038917899741

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

