



Product description

Order number	95.956
Product description	Felt marker, red, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

Product characteristics

Colour of print/label	white
Example of use	for labelling plastic tubes and reaction tubes

Size

Diameter	10 mm
Length of product	142 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Felt marker, red, waterproof

Page 2

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	10
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	300 mm
Width of case	291 mm
Height of case	392 mm
Case volume	0.0342 cbm
Weight of product	0.009 kg
Weight of case	10 kg
EAN of inner box	4038917084659
EAN case	4038917003421

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT