



Product description

Order number	70.1183.001
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 250 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	Gilson
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	2500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2500
Piece(s) / pallet	25000
Depth of case	298 mm
Width of case	464 mm
Height of case	698 mm
Case volume	0.0966 cbm
Weight of product	0.0026 kg
Weight of case	7.3 kg
EAN of inner box	4038917126441
EAN case	4038917574587

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

