



### **Product description**

Order number	70.1183.002
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 250 piece(s)/bag

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	2500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2500
Piece(s) / pallet	25000
Depth of case	298 mm
Width of case	464 mm
Height of case	698 mm
Case volume	0.0966 cbm
Weight of product	0.0026 kg
Weight of case	7.3 kg
EAN of inner box	4038917126458
EAN case	4038917574594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

