



### **Product description**

Order number	70.1186.100
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Max. volume	1250 µl
Length of product	103 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	3840
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	3840
Piece(s) / pallet	46080
Depth of case	264 mm
Width of case	475 mm
Height of case	500 mm
Case volume	0.0627 cbm
Weight of product	0.0022 kg
Weight of case	9.05 kg
EAN of inner box	4038917126489
EAN case	4038917080347

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

