




Product Specification
**Pipette tip, 1,250 µl, transparent, Biosphere® , 96
piece(s)/box**

Page 1



Product description

Order number	70.1186.200
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® , suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Pipette tip, 1,250 µl, transparent, Biosphere® , 96
piece(s)/box**

Page 2

Size

Max. volume	1250 µl
Length of product	103 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	Biosphere®
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	768
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	768
Piece(s) / pallet	79872
Depth of case	270 mm
Width of case	121 mm
Height of case	413 mm
Case volume	0.0135 cbm
Weight of product	0.0026 kg
Weight of case	2 kg
EAN of inner box	4038917160377
EAN case	4038917160360

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

