



Product description

Order number	70.1186
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume	1250 µl
Length of product	103 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	3840
Type of smallest subpackaging	bag
Piece(s) / inner box	384
Piece (s) / outer case	3840
Piece(s) / pallet	134400
Depth of case	347 mm
Width of case	258 mm
Height of case	443 mm
Case volume	0.0397 cbm
Weight of product	0.0011 kg
Weight of case	4.57 kg
EAN of inner box	4038917225236
EAN case	4038917225229

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

