



Product description

Order number	70.1187.001
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 100 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	160.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	20000
Depth of case	467 mm
Width of case	300 mm
Height of case	400 mm
Case volume	0.056 cbm
Weight of product	0.0038 kg
Weight of case	6.6 kg
EAN of inner box	4038917574617
EAN case	4038917574600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

