



Product description

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 70.3010.210 |
| Product description | Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box |

Product characteristics

| | |
|------------------------|-------------------------------------------------------------------------------------------------|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design |
| Volume of work | 10 µl |
| Filter | no |
| Coding colour | grey |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 10 µl, transparent, PCR Performance
Tested, Low retention, 96 piece(s)/box

Page 2

Size

| | |
|-------------------|-------|
| Max. volume | 10 µl |
| Length of product | 31 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | box |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 96000 |
| Depth of box | 330 mm |
| Width of box | 111 mm |
| Height of box | 148 mm |
| Depth of case | 348 mm |
| Width of case | 325 mm |
| Height of case | 245 mm |
| Case volume | 0.0277 cbm |
| Weight of product | 0.0009 kg |
| Weight of case | 5.16 kg |
| EAN of inner box | 4038917894753 |
| EAN case | 4038917894760 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT