




**Product Specification**  
**Filter tip, 10 µl, transparent, Biosphere® plus,**  
**Low retention, 96 piece(s)/box**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 70.3010.275   |
| Product description | Filter tip, working volume: 10 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design, 96 piece(s)/box |

**Product characteristics**

|                |  |
|----------------|--|
| Graduation     | yes  |
| Suitable for   | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design |
| Volume of work | 10 µl  |
| Filter         | yes  |
| Coding colour  | grey   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, 10 µl, transparent, Biosphere® plus,**  
**Low retention, 96 piece(s)/box**

**Page 2**

**Size**

|                   |       |
|-------------------|-------|
| Max. volume       | 10 µl |
| Length of product | 31 mm |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Filter material   | Polyethylene (PE)  |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | Biosphere® plus               |
| Sterilisation    | Gas sterilisation             |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | box           |
| Piece(s) / inner box          | 480           |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 96000         |
| Depth of box                  | 330 mm        |
| Width of box                  | 111 mm        |
| Height of box                 | 148 mm        |
| Depth of case                 | 348 mm        |
| Width of case                 | 325 mm        |
| Height of case                | 245 mm        |
| Case volume                   | 0.0277 cbm    |
| Weight of product             | 0.0009 kg     |
| Weight of case                | 5.16 kg       |
| EAN of inner box              | 4038917894623 |
| EAN case                      | 4038917894630 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**