



### **Product description**

Order number	70.3010.355
Product description	Filter tip, working volume: 10 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

### **Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	10 µl
Filter	yes
Coding colour	grey

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Max. volume	10 µl
Length of product	31 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

## Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	115200
Depth of box	291 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	437 mm
Width of case	307 mm
Height of case	170 mm
Case volume	0.0228 cbm
Weight of product	0.0006 kg
Weight of case	5.16 kg
EAN of inner box	4038917894838
EAN case	4038917894845

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

