



### Product description

Order number	70.3020.205
Product description	Pipette tip, working volume: 20 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	20 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 20 µl, transparent, Biosphere® plus,**  
**96 piece(s)/box**

**Page 2**

**Size**

Max. volume	20 µl
Length of product	46 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	96000
Depth of box	330 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0277 cbm
Weight of product	0.0009 kg
Weight of case	5.16 kg
EAN of inner box	4038917894227
EAN case	4038917894234

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**